



## **SDS-PAGE Protein Loading Buffer 2X (Reducing)**

**Catalog number: AR0131-20**

Boster's SDS PAGE Sample Buffer 2X (Reducing) is the most commonly used sample buffer for Sodium Dodecyl Sulfate - Polyacrylamide Gel Electrophoresis (SDS-PAGE) of denatured proteins in the Laemmli SDS-PAGE system.

This package insert must be read in its entirety before using this product. For research use only. Not for use in diagnostic procedures.

## SDS-PAGE Protein Loading Buffer 2X (Reducing)

Catalog Number: AR0131-20

### Overview

|                           |   |
|---------------------------|---|
| <b>Product Name</b>       | SDS-PAGE Protein Loading Buffer 2X (Reducing)   |
| <b>SKU/Catalog Number</b> | AR0131-20   |
| <b>Form</b>               | Liquid  |
| <b>Size</b>               | 20mL  |
| <b>Contents</b>           | 4% SDS, 20% glycerol, 200mM DTT, 0.01% bromphenol blue and 0.1 M Tris HCl, pH 6.8   |
| <b>Description</b>        | SDS PAGE Sample Buffer 2X (Reducing) is the most commonly used sample buffer for Sodium Dodecyl Sulfate - Polyacrylamide Gel Electrophoresis (SDS-PAGE) of denatured proteins in the Laemmli SDS-PAGE system. |
| <b>Equivalent</b>         | Millipore Sigma (Product No. S3401); Bio-Rad (Product No. 1610737)  |
| <b>Cite This Product</b>  | SDS-PAGE Protein Loading Buffer 2X (Reducing) (Boster Biological Technology, Pleasanton CA, USA, Catalog # AR0131-20)   |
| <b>Storage</b>            | Upon receipt store at -20°C. It is stable at -20°C for one year. Product is shipped on ice.   |

### Assay Principle

Using bromophenol blue dye, SDS-PAGE Protein Loading Buffer is a ready-to-use 2X solution. It contains 4% SDS, 20% glycerol, 200mM DTT, 0.01% bromphenol blue and 0.1 M Tris HCl. It can be used for SDS-PAGE protein loading of conventional proteins. It is especially formulated for protein sample preparation to be used in the Laemmli SDS-PAGE system. It is sufficient to prepare 12 mL protein samples. And the whole sample preparation process can be finished within 5 minutes.

## Assay Protocol

---

1. Equilibrate SDS-PAGE Protein Loading Buffer 2X to room temperature or thaw Loading Buffer in a water bath no higher than 30°C.
2. Mix one volume of SDS-PAGE Protein Loading Buffer 2X with one volume of protein sample (i.e. add 1mL protein sample into 1 mL Loading Buffer).
3. Boil sample for 3-5 min.  
**Note:** If semitransparent viscous substance remains after boiling, boil sample for another 5-10 min or add 1X SDS-PAGE Protein Loading Buffer which is diluted from 2X SDS-PAGE Protein Loading Buffer to the sample and then boil for another 3-5 min.
4. Allow sample to cool to room temperature.
5. Load sample into the wells of the SDS-PAGE gel and begin electrophoresis.
6. Stop electrophoresis when bromophenol blue dye front reaches to the bottom of the gel.

## General Notes

---

1. SDS-PAGE Protein Loading Buffer 2X contains DTT which has a slightly irritating odor.
2. SDS-PAGE Protein Loading Buffer 2X should be completely dissolved before use.
3. Please wear the lab coat and disposable gloves and operate in a laboratory hood.