

AffiniPure Goat Anti-Mouse IgG (H+L)

Basic Information

Catalog Number	BA1038-1
Size	1 mg
Ig Type	IgG
Immunogen	Mouse IgG (whole molecular)
Host	Goat
Purity	The antibody was purified from antiserum by immunoaffinity chromatography.
Form	Lyophilized
Contents	1 mg of affinity purified specific antibody (purification>99%), 0.01M PBS
Specificity	This antibody is specific for mouse IgG and shows no cross-reactivity with human/bovine/rabbit IgG.

Usage Method

Recommended Dilutions: Optimal working dilutions must be determined by end user. Working solution can not be preserved for long after being diluted with 0.01M PBS buffer solution (pH7.0-7.6) or 1% sodium acetate unless adding antiseptic as 0.05% NaN₃ and BSA to 1%.

When labeling the antibody with HRP or FITC or BIOTIN or others, reagents must be without BSA and NaN₃.

Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing.