

## AffiniPure Goat IgG (H+L)

### Basic Information

<b>Catalog Number</b>	BA1044-1
<b>Size</b>	1 mg
<b>Ig Type</b>	IgG
<b>Immunogen</b>	Goat IgG (whole molecular)
<b>Host</b>	Goat serum
<b>Purity</b>	The antibody was purified from antiserum by Protein G
<b>Form</b>	Lyophilized
<b>Contents</b>	1 mg of affinity purified specific antibody (purification>99%), 0.01M PBS

### Usage Method

Recommended Dilutions: Optimal working dilutions must be determined by end user. Working solution can not be preserved for long after being diluted with 0.01M PBS buffer solution (pH7.0-7.6) or 1% sodium acetate unless adding antiseptic as 0.05% NaN<sub>3</sub> and BSA to 1%.

When labeling the antibody with HRP or FITC or BIOTIN or others, reagents must be without BSA and NaN<sub>3</sub>.

### Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing.