

DyLight 550 Conjugated Avidin secondary antibody

Basic Information

Catalog Number	BA1134
Size	0.5 mL
Conjugation	DyLight 550
Fluorophore	Amax=562nm; Emax=576nm
Form	Concentrated, Liquid
Tested Application	IF, Flow Cytometry
Contents	0.5 mg of DyLight 550 conjugated Avidin; 0.01M PBS, 5 mg/ml BSA, 50% glycerol.
Specificity	This DyLight 550 conjugated antibody is specific for Avidin.

Suggested Working Concentration

Immunofluorescence: 1.25-5 µg/mL

Flow Cytometry: 1.25-5 µg/mL

Optimal working dilutions must be determined by end users.

Application Notes

Please centrifuge before use and dilute with neutral PBS or TBS.

The diluted solution should be used on the same day. Any unused solution should be discarded.

Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.