

## Certificate of Analysis

Product Name: Polyclonal Anti-Danio rerio (Zebrafish) kif5aa frameshift Antibody

Catalog Number: DZ41430

### Description

<b>Lot No.</b>	23DZ559G25
<b>Size</b>	200µl/vial
<b>Form</b>	Liquid
<b>Concentration</b>	500µg/ml
<b>Ig type</b>	Rabbit IgG
<b>Purification</b>	Immunogen affinity purified.
<b>Contents</b>	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN <sub>3</sub> .

### Data of Indirect ELISA

Coating concentration of antigen	0.1µg/ml	0.1µg/ml	0.1µg/ml	0.1µg/ml
Working concentration of anti-kif5aa frameshift antibody	1µg/ml	0.1µg/ml	0.01µg/ml	0.001µg/ml
Working concentration of HRP conjugated goat anti-rabbit IgG (Catalog # BA1054)	0.2µg/ml	0.2µg/ml	0.2µg/ml	0.2µg/ml
O.D. values(450nm)	2.422	2.375	0.888	0.152
Antibody Status	Pass			


### Data Evaluation Criteria

Coating concentration of antigen	0.1µg/ml	0.1µg/ml
Working concentration of the primary antibody	1µg/ml	1µg/ml
Working concentration of HRP conjugated goat anti-rabbit IgG	0.2µg/ml	0.2µg/ml
O.D. values (450nm)	≥ 1.0	<1.0
Antibody's affinity for the antigen	Strong	Weak
Antibody Status	Pass	Fail

### Storage

At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.

Approved by:



Ms. Yangjing Yue

Chief technician

Quality Control System of ELISA

Date: July 26, 2023