

Trusted by **10,000+** Scientists over 20 years




25 LONG YEARS AND COUNTING


We have spent the last 2 decades perfecting our techniques and technology to produce the highest quality biological research reagents.

BOSTER TROUBLESHOOTING EBOOK SERIES


www.bosterbio.com/ebooks



16,000+
ANTIBODIES




PICOKINE® ELISA KITS
1000+ AVAILABLE



RABBIT MONOCLONAL
ANTIBODIES 1300+



20,000+
PUBLICATIONS



Source: Biocompare.com

"This Antibody Works Perfectly!" ★★★★★ *Excellent*

SUMMARY: **Anti-GFAP Antibody (PA1239)**

"I used it for IHC on frozen sections at a dilution of 1:500. It did not need several trials to optimize the protocol. No bad things overall. I will purchase it again."

-- *Maria Teresa Dell'Anno, NEUROSCIENCE, YALE UNIVERSITY, POSTDOC*

BOSTER
antibody and ELISA experts

For more details visit, www.bosterbio.com

3942 B Valley Ave, Pleasanton, CA, 94566

Toll free: 888-466-3604 | Email: support@bosterbio.com